






### Without needles

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
	2	3/0	13 x 60 cm	black		SI30WN60-13
	3	2/0	13 x 60 cm	black		SI20WN60-13
	3.5	0	13 x 60 cm	black		SI00WN60-13
	4	1	13 x 60 cm	black		SI01WN60-13
	2	3/0	12 x 75 cm	black		SI30WN75-12
	3	2/0	12 x 75 cm	black		SI20WN75-12






### Round bodied needle – taper point – 1/2 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HR16 16 mm 	1.5	4/0	75 cm	black		SI40HR16
HR20 20 mm 	1	5/0	75 cm	black		SI50HR20
	2	3/0	75 cm	black		SI30HR20
	3	2/0	75 cm	black		SI20HR20
HR26 26 mm 	2	3/0	75 cm	black		SI30HR26
	3	2/0	75 cm	black		SI20HR26
HR30 30 mm 	3.5	0	75 cm	black		SI00HR30
	4	1	75 cm	black		SI01HR30
HR37 37 mm 	3.5	0	75 cm	black		SI00HR37



### Reverse cutting needle – 3/8 circle

NEEDLE	EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS16 16 mm 	0.7	6/0	75 cm	black		SI60DS16
	1	5/0	75 cm	black		SI50DS16
	1.5	4/0	75 cm	black		SI40DS16
	2	3/0	75 cm	black		SI30DS16
	3	2/0	75 cm	black		SI20DS16
DS19 19 mm 	1	5/0	75 cm	black		SI50DS19
	1.5	4/0	75 cm	black		SI40DS19
	2	3/0	75 cm	black		SI30DS19
	3	2/0	75 cm	black		SI20DS19
DS24 24 mm 	2	3/0	75 cm	black		SI30DS24
	3	2/0	75 cm	black		SI20DS24
	3.5	0	75 cm	black		SI00DS24



**Reverse cutting needle – 3/8 circle**

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DS30 30 mm		1.5	4/0	75 cm	black		SI40DS30
		2	3/0	75 cm	black		SI30DS30
		3	2/0	75 cm	black		SI20DS30
		3.5	0	75 cm	black		SI00DS30
		4	1	75 cm	black		SI01DS30
DS37 37 mm		3.5	0	75 cm	black		SI00DS37
		4	1	75 cm	black		SI01DS37
		5	2	75 cm	black		SI02DS37
DS45 45 mm		3.5	0	75 cm	black		SI00DS45
		4	1	75 cm	black		SI01DS45
		5	2	75 cm	black		SI02DS45
DS60 60 mm		4	1	75 cm	black		SI01DS60



**Reverse cutting needle – 1/2 circle**

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HS17 17 mm		1.5	4/0	75 cm	black		SI40HS17
		2	3/0	75 cm	black		SI30HS17
HS20 20 mm		1.5	4/0	75 cm	black		SI40HS20
		2	3/0	75 cm	black		SI30HS20
HS22 22 mm		1	5/0	75 cm	black		SI50HS22
		1.5	4/0	75 cm	black		SI40HS22
		2	3/0	75 cm	black		SI30HS22
		3	2/0	75 cm	black		SI20HS22
HS26 26 mm		1.5	4/0	75 cm	black		SI40HS26
		2	3/0	75 cm	black		SI30HS26
		3	2/0	75 cm	black		SI20HS26



**Round bodied needle – tapercut point – 1/2 circle**

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
HRT18 18 mm		1.5	4/0	75 cm	black		SI40HRT18
		2	3/0	75 cm	black		SI30HRT18



**Round bodied needle – tapercut point – 3/8 circle**

NEEDLE		EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER
DRT18 18 mm		1.5	4/0	75 cm	black		SI40DRT18
		2	3/0	75 cm	black		SI30DRT18